

# CleanSpray® HandShower+Filter AquaMedix Part Number 1415 (White) & 1302 (Chrome)

**Protect** your institution from *Legionella* and other potentially fatal waterborne bacterial infections. AquaMedix filters provide reliable bacteria and protozoan cyst reduction for safer and cleaner bathing and shower water for visitors, patients, and employees without incorporating chemicals into the system.

## CleanSpray® HandShower Contains Hollow Fiber Micro-Filter 0.15 µm Membrane:

- ✓ Absolute 0.15 micron particle retention for bacteria reduction such as *Legionella* and *Pseudomonas*
- ✓ 99.99999% bacteria reduction (> 7 Log)
- ✓ Captures contaminants in water before they are aerosolized
- ✓ FDA - FDA Class I medical device
- ✓ EPA – EPA Pest control device
- ✓ Can be used in combination with secondary water disinfection chemical treatment systems
- ✓ Environmentally friendly alternative handshower with replaceable filter
- ✓ 3-month (#1413) or 1-month (#1430) replacement filters



**For use in:** Hospitals, Medical Facilities, Adult Care Facilities, Health Clubs, Homes, Swimming Pools, etc.

Technical data	CleanSpray® Handshower
Membrane, hollow fibers microfilter	Absolute 0.15 micron
Bacteria retention*	99.99999% (>7 log)
Maximum life†	3-month†
Estimated capacity volume†	525 gal† to 1800 gal†
Initial flowrate† (at 72 psi)	5.80 gpm
Maximum operating pressure	Pressure 25-88 psi, Max. 88 psi
Operating temperature	40-110 °F, Max. temp 140 °F
Membrane material	PES, (Polyethersulfone) proprietary
Cartridge dimensions	1-3/16 in. X 6 in.
Housing dimensions	3.50 in. X 10.25 in. high
Part number	1415 (white) & 1302 (chrome)
Connection	Standard ½" NPT



CleanSpray® HandShower+Filter  
Part #1415 White



CleanSpray® HandShower+Filter  
Part #1302 Chrome

\*Bacteria retention tested to 7 log *Br. Dimunta*

†Assuming a potable water supply quality in accordance with the Safe Drinking Water Act  
FDA Class I Medical Device/EPA Pest Control Device | EPA Est. No. 89767-MN-001 | U.S. Patent No. 5,240,862

AquaMedix, LLC | 7169 Shady Oak Road, Minneapolis, MN 55344 | [www.aquamedix.net](http://www.aquamedix.net) | [info@aquamedix.net](mailto:info@aquamedix.net) | 952-479-0636

Manufactured and Distributed in the United States by AquaMedix LLC.



Filters are intended for use on visually clear drinking water that meets all other public health standards. Filters are not intended for reducing pathogenic (disease-causing) virus particles. The Micro-Filter is a USEPA Pest Control Device - EPA Est. No. 89767-MN-001, U.S. Patent No. 5,240,862. **Potable water service only.** Filter not for use with any injection or infusion applications or applications requiring United States Pharmacopeia sterile water. CleanSpray® is a Registered Trademark of AquaMedix, LLC, Reg. No. 4,450,320. All rights reserved. | CleanSpray\_Handshower\_HSw&C\_rev102518 Legionella filters Made in Belgium