

Filtration technologies for immediate control of *Legionella* and other waterborne bacteria

AquaMedix Point-of-Use Products



Part #	Product Name	Filter max. life†	Micron	FDA/EPA Registration	
1690	CleanSpray Faucet Filter FX (3M)	3 months	0.2	FDA/EPA	1 or 3 month 0.20 micron filters. 55/64"-27 F connector enables fast and easy replacements and obscures the unlock mechanism to prevent removal by users minimizing the possibility of contamination.
1630	CleanSpray Faucet Filter FX (1M)	1 month	0.2	FDA/EPA	
Each Faucet Filter requires a Compliance Connector to connect to faucet.					
48116	CleanSpray Faucet Compliance Connector 55/64"-27 F	Compliance Connector includes a 15/16"-27 M x 55/64"-27 M adapter and extra flat gasket seal.			
1413	HandShower Replacement Filters (3M)	3 months	0.2	FDA/EPA	1 or 3 month 0.20 micron filters mount inside chrome or white shower handles. Handles equipped with check valve and flow controller to meet DOE flow requirements. Shower handles have a standard 1/2" male pipe thread inlet connections.
1430	HandShower Replacement Filter (1M)	1 month	0.2	FDA/EPA	
1302	CleanSpray HSC HandShower+Filter, Chrome	3 months	0.2	FDA/EPA	
1415	CleanSpray HSW HandShower+Filter, White	3 months	0.2	FDA/EPA	
1090W	Shower Pre-Filter System, White	6 month	2.0		Pre-filter removes larger particles in the water preventing filter inside the shower handle from prematurely clogging.
1091	Pre-Filter Replacement Cartridge	6 month	2.0		
For use in standard 10" inch housings				EPA	For use in standard 10" housings. 2.5 GPM at 30 PSI (max. flowrate 10 GPM). Applications: Inline, Ice Machines, etc. Micro 20" limited available.
1042	AquaMedix Micro 10 (inch)	6 months	0.2	EPA	
1379	AquaMedix Micro 20 (inch) 20" Special order	6 months	0.2	EPA	
1009	AquaMedix Dual Stage CL CarbonMicro	6 months	0.2	EPA	Carbon shell pre-filter (5 µm) eliminates chlorine, heavy metals, drug residues bad taste and other contaminants. Micro-Filter removes bacteria.



Baclyser® by Aqua-free

Part #	Product Name	Filter max. life†	Micron	FDA/EPA Registration	
49010	Baclyser TL (2M) Tap	2 months	0.2	FDA/EPA	The disposable TL with laminar outlet and TR with shower rose outlet are individually barcoded for traceability. Filters are easily changed via the quick release coupling.
49110	Baclyser TL (1M) Tap	1 month	0.2	FDA/EPA	
49011	Baclyser TR (2M) Tap or Shower	2 months	0.2	FDA/EPA	
49111	Baclyser TR (1M) Tap or Shower	1 month	0.2	FDA/EPA	
Quick-release coupling adapter required					
49012	Baclyser S (2M) Shower	2 months	0.2	FDA/EPA	The disposable S shower has a copper-plated shower rose to help prevent retrograde contamination. Each individually barcoded for traceability.
49112	Baclyser S (1M) Shower Connection 1/2" Male Pipe	1 month	0.2	FDA/EPA	
49414	AS WallShower Connection 1/2" Male Pipe	4 months	0.2	FDA/EPA	The disposable AS WallShower is for use with shower arms that come directly out of the wall.
Inline filters					
49404	Baclyser IL3 Filter	6 months	0.2	FDA/EPA	High flow rate, high capacity multi-use for use with ice machines, fixed showers, under sink, drinking fountain etc. Quick Connect fittings for easy change outs available.
49406	Baclyser IL5 Filter	6 months	0.2	FDA/EPA	
49408	Baclyser IL10 Filter Connection 1/2" Male Pipe. See connector kits.	6 months	0.2	FDA/EPA	

Contact your distributor or AquaMedix, LLC for product pricing and volume discounts. A complete list of all AquaMedix filtration products and connection parts available upon request. Installation kits and instruction guides are available for unique and special plumbing connections.

AquaMedix, LLC | 7169 Shady Oak Rd. Eden Prairie, MN 55344 | 952-479-0636

www.aquamedix.net | info@aquamedix.net

CleanSpray and Baclyser filters are FDA Class I Medical Devices/ EPA Pest Control Devices. EPA Est. No. 89767-MN-001

†Assuming a water supply in accordance with the Safe Drinking Water Act (as a guide value).

For terms and conditions see our website

Updated 05/15/20

AquaMedix Filtration Product list rev051520 7.94

	Part #	Product Name	
Installation kits with quick-connect and other connector parts	QCK01	Baclyser InLine Comprehensive Installation Kit	
	QCK03	Baclyser InLine Series-40 New Installation Kit	AquaMedix Quick-Connect Kits are designed to simplify and speed the installation and future change out of Baclyser Inline filters. Kits are available for a variety of plumbing options to identify the best approach for new installations, retrofits, or compatibility with other manufacturer specialized install kits.
	QCK04	Baclyser InLine Series-40 Replacement Kit	
	QCK05	Baclyser InLine Series-50 Replacement Kit	
	QCK06	Baclyser InLine Series-50 New Installation Kit	
	QCK09	Baclyser Inline Screw-On Push Connect Kit	
	QCK16	Baclyser Faucet Tap Connect Kit	Faucet tap connection kits provide a selection of the most commonly used adapter parts necessary for connecting to a faucet.
QCK17	CleanSpray Faucet Tap Connect Kit		

	Part #	Description
	48112*	Quick-Coupling, 55/64"-27 Female
	48114*	Quick-Coupling, 15/16"-27 Male
	CRC050	1/2" F x 1/2" F, IPS Copper Plated Brass (CPB) Coupling (Lead Free) Alt. part # 341404
	15 3360 5	NeoPerl 15/16"-27 Male x 55/64"-27 Male
	15 3410 5	NeoPerl 3/4"-27 Female x 15/16"-27 Male
	15 3390 5	NeoPerl 13/16"-27 Male x 55/64"-27 Male
	15 3110 5	NeoPerl 1/2"-14 FPT x 55/64"-27 Male
	40CB-S1-06	Quick-Coupling, 3/8" MPT
	40ACV-SB1-06	Quick-Coupling, Valved 3/8" MPT, plastic
	40AC-PB1-06	Quick-Connect plug 3/8" MPT, plastic
	PLC29006	CPC Quick-Disconnect Coupling Insert, plug x 3/8" OD; Non-Valved
	BBFM-1238	1/2" MPT x 3/8" FPT Lead Free Brass Bushing
	AFAB0607CW	DMFit Female Adapter 3/8" Tube x 1/2" BSPP thread – Cone Seal Alt. part # PI451214FS*
	AFAB0407CW	DMFit Female Adapter - BSP thread 1/4" Tube x 1/2" BSPP Thread – Cone Seal
	PSEI6012U9	John Guest 3/8" Tube x 3/8" Compression (9/16"-24 UNEF)
	PI221212	John Guest Stem Elbow Connector - 3/8" Stem x 3/8" Tube
	PI4512F6S	John Guest Female Adapter NPTF - 3/8" x 3/8" NPTF
	C-2500W	2.5" Single Filter Mounting Bracket, White Alt. part # PPC202W
	MS1010	#10-3/4" Screw (for C-2500W Mounting Bracket)
	DPE 06-B	3/8" OD Tubing, NSF (Sold per foot)
	430-073*	Thread coupling 1/2" x 3/8" FxF NPT, plastic* Alternative part # L805100*
	430-005	Thread coupling 1/2" FxF NPT, plastic Alternative part # L805000

*Available through AquaMedix, or special order with other suppliers.

Contact your distributor or AquaMedix, LLC for product pricing and volume discounts. A complete list of all AquaMedix filtration products and connection parts available upon request. Installation kits and instruction guides are available for unique and special plumbing connections.

AquaMedix, LLC | 7169 Shady Oak Rd. Eden Prairie, MN 55344 | 952-479-0636

www.aquamedix.net | info@aquamedix.net

CleanSpray and Baclyser filters are FDA Class I Medical Devices/ EPA Pest Control Devices. EPA Est. No. 89767-MN-001

Updated 05/15/20

For terms and conditions see our website

AquaMedix Filtration Product list rev051520 7.94