







# Filtration technologies for immediate control of *Legionella* and other waterborne bacteria





## AquaMedix Point-of-Use Products



|  | Part #                                       | Product Name                                       | Filter max. life†   | Micron | FDA/EPA Registration |   |
|--|--|--|---|--------|----------------------|---|
|  | 1690   | CleanSpray Faucet Filter 3 month                   | 90 day  | 0.20   | Pending              | 1 or 3 month 0.20 micron filters. 55/64"-27 F connector enables fast and easy replacements and obscures the unlock mechanism to prevent removal by users minimizing the possibility of contamination.                         |
|  | 1630   | CleanSpray Faucet Filter 1 month                   | 31 day  | 0.20   | Pending              |   |
|  | 48116  | CleanSpray Faucet Compliance Connector 55/64"-27 F | Compliance Connector includes a 15/16"-27 M x 55/64"-27 M adapter and extra flat gasket seal. |        |                      |   |
|  | 1413   | HandShower Replacement Filters (3M)                | 90 day  | 0.20   | FDA/EPA              | 1 or 2 month 0.20 micron filters mount inside chrome or white shower handles. Handles equipped with check valve and flow controller to meet DOE flow requirements. Shower handles have a standard 1/2" male pipe thread inlet |
|  | 1430   | HandShower Replacement Filter (1M)                 | 31 day  | 0.20   | FDA/EPA              |   |
|  | 1302   | CleanSpray HandShower+Filter, Chrome               | 90 day  | 0.20   | FDA/EPA              |   |
|  | 1415   | CleanSpray HandShower+Filter, White                | 90 day  | 0.20   | FDA/EPA              |   |
|  | 1090W  | Shower Pre-Filter System, White                    | 6 month   | 2.0    |                      | Pre-filter removes larger particles in the water preventing filter inside the shower handle from prematurely clogging.  |
|  | 1091   | Pre-Filter Replacement Cartridge                   | 6 month   | 2.0    |                      |   |
|  | <b>For use in standard 10" inch housings</b> |  |   |        | EPA                  | For use in standard 10" housings. 2.5 GPM at 30 PSI (max. flowrate 10 GPM). Applications: Inline, Ice Machines, etc. Micro 20" available as special order only  |
|  | 1042   | AquaMedix Micro 10 (inch)                          | 6 month   | 0.20   | EPA                  |   |
|  | 1379   | AquaMedix Micro 20 (inch) 20" Special order        | 6 month   | 0.20   |                      |   |
|  | 1009   | AquaMedix Dual Stage CL CarbonMicro                | 6 month   | 0.20   | EPA                  | Carbon shell pre-filter (5 µm) eliminates chlorine, heavy metals, drug residues bad taste and other contaminants. Micro-Filter removes bacteria.  |



### Baclyser® by Aqua-free

|  |       | Product Name                                   | Filter max. life† | Micron | FDA/EPA Registration |   |
|--|-------|--|-------------------|--------|----------------------|---|
|  | 49010 | Baclyser TL (2M) Tap                           | 62 day            | 0.20   | FDA/EPA              | The disposable TL with laminar outlet and TR with shower rose outlet are individually barcoded for traceability. Filters are easily changed via the quick release coupling. |
|  | 49110 | Baclyser TL (1M) Tap                           | 31 day            | 0.20   | FDA/EPA              |   |
|  | 49011 | Baclyser TR (2M) Tap or Shower                 | 62 day            | 0.20   | FDA/EPA              |   |
|  | 49111 | Baclyser TR (1M) Tap or Shower                 | 31 day            | 0.20   | FDA/EPA              |   |
|  |       | Quick-release coupling adapter required        |                   |        |                      |   |
|  | 49012 | Baclyser S (2M) Shower                         | 62 day            | 0.20   | FDA/EPA              | The disposable S shower has a copper-plated shower rose to help prevent retrograde contamination. Each individually barcoded for traceability.                              |
|  | 49112 | Baclyser S (1M) Shower                         | 31 day            | 0.20   | FDA/EPA              |   |
|  | 49414 | AS WallShower                                  | 4 month           | 0.20   | FDA/EPA              | The disposable AS WallShower is for use with shower arms that come directly out of the wall.  |
|  |       | <b>Inline filters</b>                          |                   |        |                      |   |
|  | 49404 | Baclyser IL3 Filter                            | 6 month           | 0.20   | FDA/EPA              | High flow rate, high capacity multi-use for use with ice machines, fixed showers, under sink, drinking fountain etc. Quick Connect fittings for easy change outs available. |
|  | 49406 | Baclyser IL5 Filter                            | 6 month           | 0.20   | FDA/EPA              |   |
|  | 49408 | Baclyser IL10 Filter                           | 6 month           | 0.20   | FDA/EPA              |   |
|  |       | Connection 1/2" Male Pipe. See connector kits. |                   |        |                      |   |

Contact your distributor or AquaMedix, LLC for product pricing and volume discounts. A complete list of all AquaMedix filtration products and connection parts available upon request. Installation kits and instruction guides are available for unique and special plumbing connections.

AquaMedix, LLC | 7169 Shady Oak Rd. Eden Prairie, MN 55344 | 952-479-0636

[www.aquamedix.net](http://www.aquamedix.net) | [info@aquamedix.net](mailto:info@aquamedix.net)

AquaMedix Product list rev091019

CleanSpray and Baclyser filters are FDA Class I Medical Devices/ EPA Pest Control Devices. EPA Est. No. 89767-MN-001

†Assuming a water supply in accordance with the Safe Drinking Water Act (as a guide value).